



Answers for Science.
Knowledge for Life.™




MarkerView™ Software 1.3

New Features

For Research Use Only. Not for use in diagnostic procedures. RUO-MKT-11-4893-A

MarkerView™ Software Overview

- A data visualization tool designed for scientists who wish to visualize their data in terms of sample groupings and apply statistics in order to gain valuable insight into any trends within their mass spectral data.
- MarkerView is unique in that SCIEX users can explore statistical correlations with direct connections back to the raw data. This allows them to find meaningful relationships much more quickly.
- Target Applications:
 - Metabolomics
 - Lipidomics
 - Proteomics
 - Food Authenticity
 - Water Testing



MarkerView™ Software 1.3



New Features



- Support for *.wiff2 format
 - X500R QTOF System
- Import wizard improves ease of learning
- Changes to t-test view
- Box and whiskers plots
- Infusion MS/MSALL Support
- SWATH® Acquisition support
- Most likely ratio (MLR) normalisation
- Custom sample columns
- ‘Set Names’ script
- Speed and other small improvements


Support for *.wiff2 Format

- Supports *.wiff2 format enabling Metabolomics and Food/Environmental workflows on the X500R QTOF System




Import Wizard

- A wizard enables several steps to be combined:
 - Selecting input files
 - Peak finding and alignment parameters
 - Sample group definition (previously set in Samples Table)
 - Point symbols for groups (previously set in Options)
 - Automatic processing such as PCA or t-test (previously done explicitly)



Updated T-tests



Unequal variance test - when the two groups have different variances and/or group sizes

Sorts automatically by p-value



1.3 - [Box and Whiskers Plot]

View Analyze Window Help

Compare: 1 to blank n1 = 3, n2 = 9

Row	Index	Peak Name	m/z	Ret. Time	Group	Use	t-value	p-value	Mean
1	106	323.2/13.0 (106)	323.1720	12.95	1	<input checked="" type="checkbox"/>	29.48	4.7085e-11	2.966e3
2	70	266.1/12.8 (70)	266.1189	12.77	1	<input checked="" type="checkbox"/>	28.31	7.0192e-11	1.689e3
3	63	252.1/12.8 (63)	252.1008	12.76	1	<input checked="" type="checkbox"/>	26.22	1.5004e-10	9.159e2
4	84	280.1/12.7 (84)	280.1335	12.73	1	<input checked="" type="checkbox"/>	25.53	1.9480e-10	1.128e3
5	64	253.1/12.8 (64)	253.1076	12.76	1	<input checked="" type="checkbox"/>	25.39	2.0559e-10	8.027e2
6	82	279.2/12.8 (82)	279.1827	12.77	1	<input checked="" type="checkbox"/>	23.77	3.9440e-10	7.034e2


Profile Plot appears automatically and links to table



New Box & Whiskers plot appears automatically



Box and Whiskers Plots

- Standard way of visualizing statistical data across groups
- Spacing between the different parts of the box indicate the degree of spread and skewness in the data
- Highlights outliers
- Automatically generated after t-test



Infusion MS/MS^{ALL} Support

- Infusion MS/MS^{ALL} acquisition is a technique in which a sample is infused and MS/MS is acquired for each precursor (at unit resolution) over a wide mass range
 - Typically used for lipidomics





SWATH® Acquisition Support

- Import raw SWATH acquisition datafiles
- MS/MS spectra can be viewed for selected features (similar to IDA functionality).
- Targeted Spectrum Finder removes fragments with mismatching LC peak profiles


Peak Name	m/z	Ret. Time	Group	Use	Accum80 S1A rep4
			(Monoisotopic)	<input checked="" type="checkbox"/>	2.566e4
			(Isotope)	<input checked="" type="checkbox"/>	2.326e4
			(Monoisotopic)	<input checked="" type="checkbox"/>	2.268e4
			(Monoisotopic)	<input checked="" type="checkbox"/>	2.169e4
					160e4
					139e4
					119e4
					118e4
					114e4
					056e4
					047e4
					036e4
					887e4
					808e4

Row	Index	Peak Name	m/z	Ret. Time	Group	Use	Accum80 S1A rep
23	7592	1042.6/59.5 (7592)	1042.632	59.45	(Monoisotopic)	<input checked="" type="checkbox"/>	2.268e4
24	516	458.3/59.5 (516)	458.2575	59.54	(Monoisotopic)	<input checked="" type="checkbox"/>	2.169e4
25	3893	652.4/52.6 (3893)	652.3884	52.60	(Monoisotopic)	<input checked="" type="checkbox"/>	2.160e4
26	5046	722.0/37.5 (5046)	722.0185	37.53	(Isotope)	<input checked="" type="checkbox"/>	2.139e4
27	5421	744.0/59.5 (5421)	744.0521	59.54	(Isotope)	<input checked="" type="checkbox"/>	2.119e4




Most Likely Ratio (MLR) Normalisation

- Often need to normalize responses for a sample to allow for the fact that the *absolute* sample amount might not be constant (due to different starting amounts, sample prep differences, etc.)
- Typically used for protein/peptide normalization using large numbers of endogenous peptides (features)




For more information on the normalization strategy:

- Lambert et al. (2013) Nature Methods, 10, 1239-1245
- SCIEX Community discussion

A screenshot of a blog post from SCIEX. The title is 'What is the normalization strategy used within the Protein Expression Assembler app?'. The author is Christiel Hunter, and the post was published on February 16, 2015, at 8:13 PM. The post discusses the normalization strategy used in the Protein Expression Assembler application. It explains that normalization is utilized as the first step before computation of fold change numbers for proteins between samples. The normalization starts at the lowest level of the replicates (the biological and technical replicates together) within an experimental group. The algorithm determines the ratio for each feature between each pair of replicates within the group which creates a ratio histogram for each sample (matrix of comparisons within a single experimental group). A graph titled 'Intensity-ratio histogram for two samples' shows a bell-shaped curve representing the distribution of intensity ratios. The x-axis is labeled 'Intensity ratios (log₁₀)' and has a '0' mark. A vertical dashed line indicates the 'Most likely ratio'. A horizontal double-headed arrow indicates the 'Normalization factor'. A vertical double-headed arrow indicates the 'Inversely proportional to reproducibility'.


Custom Sample Columns

- Ability to create arbitrary columns to allow metadata entry
- Allows switching between different ways of grouping data for visualisation or supervised algorithms (t-test, PCA-DA).
- Populated manually in MarkerView™ Software or by adding custom fields to Analyst® Software batch.



'Set Names' Feature

- Utility which allows variables to be named based on their mass and retention time.
- Assuming you have a list of such masses and retention times for known compounds
 - More useful than a list of m/z-RT ion pairs
- Unknown compounds can be excluded, otherwise they remain with the default names.



Other New Features

- Peak-finding has been made 'multi-threaded'
 - On a multiple-core computer, processing will be faster (when there are multiple samples)
- 'Impact' column in PCA Loadings Table
 - Measure of the importance of each variable to separation

Trademarks/Licensing

AB Sciex is doing business as SCIEX.

For Research Use Only. Not for use in diagnostic procedures.

© 2016 AB Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.



Answers for Science.
Knowledge for Life.™

