

Identification of Ginsenosides Using the SCIEX X500R QTOF System

Wang Sha, Cheng Haiyan, Liu Ting, Li Lijun, Jin Wenhai[Author]
SCIEX, Pacific Applications Support Center (Beijing). China

Background

Ginseng is one of the most valued herbs. It has properties of nourishing vital strength, tranquilizing the mind, promoting secretions, and supplementing deficiencies. Modern medical research shows that ginseng is effective in preventing cancer, countering aging, acting as an antiarrhythmic agent, and has hypoglycemic, hypolipidemic, and immune-stimulating properties. Its main active components are ginsenosides.

Ginsenosides are triterpenoid chemical compounds; based on the glycosyl structure, they can be divided into tetracyclic triterpenes of the dammarane type and pentacyclic triterpenes of the oleanane type. The dammarane type can be divided into ginseng diols and ginseng triols. Because ginsenosides have many components, different species and sources yield differences in composition^[1], so a full identification of the ginsenoside composition and accurate analysis of its structure currently requires extensive literature and document research. At the same time, the analytical result obtained by different technologies can be quite difficult to verify with data, which can complicate quality evaluation and material basis.

The SCIEX X500R QTOF high resolution mass spectrometer requires a single injection to collect high quality MS and MS/MS data. Combined with an expansive high resolution MS/MS database of Traditional Chinese Medicine (TCM) active ingredients, SCIEX OS software automatically determines the theoretical molecular weight and isotope pattern distribution and simultaneously matches it with the MS/MS database. Comprehensive scoring allows intuitive, rapid, and accurate ginsenoside component identification.

The SCIEX high resolution database of Chinese medicine is based on “Chinese pharmacopeia” Part 1, TCM active ingredients. It includes almost a thousand compounds such as saponins, flavonoids, flavonoid glycosides, triterpenes, phenylethyl glycosides, and organic acids.

This document describes the workflow for analysis of ginsenosides in Chinese medicine using the SCIEX OS ultra-efficient data processing software and the high resolution Chinese medicine MS/MS database on the SCIEX X500R QTOF high resolution mass spectrometer system.

The software has a simple user interface, and the workflow is intuitive.

Experimental Process

1. Using TOF-IDA- mode (Top 10 MS/MS per cycle), inject a sample and simultaneously obtain primary precursor ions and secondary daughter ion information. This allows confidence in the compound identification and allows all data to be acquired in a single shot alleviating time involved in a two-injection workflow.
2. Search known ginsenoside components; according to the accurate mass, isotope distribution, and Chinese medicine matching data, identify the compounds.
3. Use Ginsenoside Rg2 as an internal standard to verify accuracy of match results.
4. Use the accurate mass, characteristic fragment ions, and relative retention times to enhance identification of ginsenoside isoforms.
5. A total of 51 commonly observed ginsenoside components have already been identified.




Figure 1 Workflow for using the SCIEX OS X500R QTOF high resolution mass spectrometer and the Chinese medicine MS/MS database to identify ginsenoside components

Sample Preparation

1. Accurately measure 5.0g ginseng powder into a 50mL centrifugation tube.
2. Add 25mL 90% methanol water, agitate 5 min.
3. Immerse in an ice bath overnight.
4. Ultrasonicate 30 min, at 4 deg. C, then centrifuge at 10000r/min for 12 min.

5. Remove the supernatant and pass through a 0.22 μ m filter.

Liquid Chromatography (LC) Conditions

Chromatographic Column: Phenomenex Kinetex C18, 2.1*100mm, 2.6 μ m;

Mobile phase: Gradient elution is used

Negative ions: A is H₂O (containing 0.05% formic acid); B is acetonitrile;

Flow rate: 0.25mL/mL

Column temperature: 40°C

Injection volume: 3 μ L

Table 1. Elution conditions

| Time (min) | A% | B% |
|------------|----|-----|
| 0 | 90 | 10 |
| 0.5 | 90 | 10 |
| 5.0 | 50 | 50 |
| 35.0 | 10 | 90 |
| 40.0 | 0 | 100 |
| 40.1 | 90 | 10 |
| 45.0 | 90 | 10 |

Mass Spectrometry Conditions

Scanning method: TOF-IDA-10 MS/MS qualitative;


ESI ion source parameters:

Air curtain gas CUR: 35psi;

IS voltage: -4500V; Source temperature: 550°C

Cone voltage: -80V;

Atomizing gas GAS1: 55psi; Auxiliary gas GAS2: 55psi



Application of SCIEX OS Software for Ginsenoside Analysis

SCIEX OS Software platform provides simultaneous mass spectrometer control, method editing, data analysis, and result reporting.

1. Data acquisition

Data acquisition is performed on extracted samples according to conditions described. The Explorer data processing options can be used to open the acquired high-resolution data, and perform any data QC as in Fig. 2.




Figure 2. Acquired high-resolution TOF MS-IDA-TOF MS/MS data. Fig. 2A shows a full TOF MS scan, and Fig. 2B is an IDA TOF MS/MS spectrogram of a sample.

2. Editing of data processing methods

Using targeted component analysis, one can enter or copy known ginsenosides to the component options, including their name and molecular composition, as shown in Fig. 3.




Figure 3. Input chemical compound list

At the same time, select the database to search (this study uses the TCM MS/MS Library) and configure the confidence levels, as shown in Fig. 4.




Figure 4. Setting SCIEX OS Software confidence levels

Confidence intervals are primarily used in compound identification and verification, including theoretical mass numbers, isotope distribution, retention times (can be omitted if unavailable), and MS/MS spectrum matching in the database. Each score is calculated based on the configuration, and an overall score is determined based on the weight of the four parameters.

3. Data processing and Results Viewing

Use built-in processing methods to processing data. Select “process,” and the software will list results based on 4 established confidence intervals. Use the “signal indicator” to easily obtain the results, including MS/MS matched results. Identification of Ginsenoside Rg2 is given as an example; results are shown in Fig.5.




Figure 5. Screening results display

Fig. A is the display of the data processing results; B is the display of an extracted ion chromatogram. Retention time is useful for compound verification. Fig. C shows the mass spectrum and isotope distributions; the upper portion is the actual measured value, the lower, gray part shows the

theoretical MS pattern, and the two match well. Fig. D is the mirror image display for the MS/MS database match results: the upper, blue portion is the MS/MS spectrum acquired for ginsenoside, and the lower, gray portion shows the MS/MS spectrum from the database. Results are compared clearly. The lower right corner shows database search results, and with an overall score above 90, matching results are excellent.

Verification of identification results

For ginsenoside in the negative ion mode, the excimer ion peaks are mainly present at $[M-H]^-$ and $[M+HCOOH-H]^-$; the structure of the saponin component is described by its secondary signature fragment loss of HCOOH and daughter ions of sugars, e.g.: -46 -162 (Loss HCOOH& Glu); -46-146 (Loss HCOOH& Rhamnose) or -46-132 (Loss HCOOH& arabinose), and fragment 161 forms readily. For ginsenoside Rg2 secondary fragment structure analysis, see Fig. 6.




Figure 6. Secondary structural analysis spectrogram for ginsenoside Rg2

Using MS/MS information, a neutral loss (NL) loss of 162 and 146, and the signature fragment for dammarane triol at m/z 475, this structure can be described as dammarane triol +glucose+rhamnose, which is similar to ginsenoside Rg2 and matches search results.

Ginsenosides are made of saponins and sugars such as glucose, rhamnose, and arabinose that may be linked at various positions on the saponin. At the same time, since saponins and sugars have different structures, they may form stereoisomers. In database searches, isomers link to secondary signature fragment ion information, retention time, and relevant literature [2]. Verification of known ginsenosides is necessary, and secondary fragment information is critical for structural identification of Chinese medicine components

Table 2. List of identified ginsenosides

| Index | Compound name | Molecular formula | [M+HCOOH-H]- | ss Error (ppm) | Retention time (min) | MS/MS |
|-------|---------------------------------|-------------------|--------------|----------------|----------------------|-----------------------------------|
| 1 | 20R-ginsenoside Rg2 | C42H72O13 | 829.4944 | 2.30 | 13.93 | m/z783, 637, 475,391,161 |
| 2 | 20S-ginsenoside Rg2 | C42H72O13 | 829.4944 | 1.78 | 14.01 | m/z783, 637, 475,391,161 |
| 3 | 20S-ginsenoside F2 | C42H72O13 | 829.4944 | 1.50 | 21.05 | m/z783, 621, 459,161 |
| 4 | 20R-ginsenoside F2 | C42H72O13 | 829.4944 | 1.10 | 20.10 | m/z783, 621, 459,161 |
| 5 | 20S-ginsenoside Rg3 | C42H72O13 | 829.4944 | -0.80 | 26.09 | m/z783, 621,459,375,161 |
| 6 | 20R-ginsenoside Rg3 | C42H72O13 | 829.4944 | -0.80 | 26.54 | m/z783, 621,459,375,161 |
| 7 | Ginsenoside Rg4 | C42H70O12 | 811.4838 | -1.20 | 7.79 | m/z765,619,457,161 |
| 8 | Ginsenoside F4 | C42H70O12 | 811.4838 | -1.70 | 12.22 | m/z765,619,457,161 |
| 9 | Ginsenoside Ic | C42H70O12 | 811.4838 | -0.90 | 28.58 | m/z765,619,457,161 |
| 10 | Ginsenoside Rk1 | C42H70O12 | 811.4838 | -1.20 | 29.04 | m/z765,603,441,161 |
| 11 | Ginsenoside Rk1 (isomer) | C42H70O12 | 811.4838 | 1.40 | 41.33 | m/z765,603,441,161 |
| 12 | Pseudoginsenoside F11 | C42H72O14 | 845.4893 | -1.60 | 12.52 | m/z799,637,475,161 |
| 13 | Ginsenoside Rf | C42H72O14 | 845.4893 | -1.90 | 8.63 | m/z799,637,475,161 |
| 14 | Ginsenoside Rf (isomer) | C42H72O14 | 845.4893 | -1.60 | 8.35 | m/z799,637,475,161 |
| 15 | Ginsenoside Rg7 | C42H72O14 | 845.4893 | -2.00 | 7.64 | m/z799,637,475,161 |
| 16 | Ginsenoside Rg1(Ginsenoside A2) | C42H72O14 | 845.4893 | -1.60 | 7.40 | m/z799,637,475,161 |
| 17 | Ginsenoside F1 | C36H62O9 | 683.4365 | 0.90 | 15.85 | m/z637,475,391,161 |
| 18 | 20(S)-Ginsenoside Rh1 | C36H62O9 | 683.4365 | 0.70 | 14.23 | m/z637,475,391,161 |
| 19 | 20(R)-Ginsenoside Rh1 | C36H62O9 | 683.4365 | 0.90 | 14.10 | m/z637,475,391,161 |
| 20 | Ginsenoside Re | C48H82O18 | 991.5472 | -0.84 | 8.47 | m/z945, 799,637,475 |
| 21 | Ginsenoside Rd | C48H82O18 | 991.5472 | -0.70 | 21.22 | m/z945,783,621,459,375,161 |
| 22 | Ginsenoside Rd (isomer) | C48H82O18 | 991.5472 | -0.70 | 22.83 | m/z945,783,621,459,161 |
| 23 | pseudo-Ginsenoside RT2 | C41H70O14 | 831.4737 | -1.50 | 7.41 | m/z785,653,491,391 |
| 24 | Ginsenoside Rb2 | C53H90O22 | 1123.5895 | -1.90 | 18.45 | m/z1077,945,783,621,459 |
| 25 | 20(S)-Ginsenoside Rc | C53H90O22 | 1123.5895 | -2.10 | 19.44 | m/z1077,945,783,621,459 |
| 26 | 20(R)-Ginsenoside Rc | C53H90O22 | 1123.5895 | -1.90 | 19.76 | m/z1077,945,783,621,459 |
| 27 | Ginsenoside Rb1 | C54H92O23 | 1153.6001 | -1.00 | 17.99 | m/z1107, 945,783,621,459 |
| 28 | 20(S)-Ginsenoside-Rh2 | C36H62O8 | 667.4416 | -0.50 | 29.52 | m/z621,459,375 |
| 29 | 20(R)-Ginsenoside-Rh2 | C36H62O8 | 667.4416 | -0.50 | 30.50 | m/z621,459,375 |
| 30 | Ginsenoside Rd+Acetylation | C50H84O19 | 1033.5578 | -1.40 | 20.26 | m/z987, 945, 928, 783, 621,459 |
| 31 | Ginsenoside Re+Acetylation | C50H84O19 | 1033.5578 | -1.52 | 20.74 | m/z987, 945, 928, 783, 621,459 |
| 32 | Pseudoginsenoside RT5 | C36H62O10 | 699.4314 | -1.10 | 7.75 | m/z699,653,491,329,161 |
| 33 | Ginsenoside Ra1 | C58H98O26 | 1255.6317 | -0.63 | 17.23 | m/z1209,1077,945,783,621,459 |
| 34 | Ginsenoside Ra2 | C58H98O26 | 1255.6317 | -0.70 | 18.40 | m/z1209,1077,945,783,621,459 |
| 35 | Chikusetsusaponin III | C47H80O17 | 961.5367 | -1.00 | 23.65 | m/z915,783,621,459,375 |
| 36 | Ginsenoside Rs2 | C55H92O23 | 1165.6001 | 1.60 | 20.17 | m/z1119,1077,1059,945,783,621,459 |
| 37 | Ginsenoside Rs2 (isomer) | C55H92O23 | 1165.6001 | 1.84 | 19.57 | m/z1119,1077,1059,945,783,621,459 |
| 38 | Ginsenoside Rs1 | C55H92O23 | 1165.6001 | 1.00 | 18.67 | m/z1119,1077,1059,945,783,621,459 |
| 39 | Ginsenoside Rs1 (isomer) | C55H92O23 | 1165.6001 | 1.00 | 17.73 | m/z1119,1077,1059,945,783,621,459 |
| 40 | Ginsenoside R1 | C47H80O18 | 977.5316 | -1.20 | 7.87 | m/z 931, 799, 637, 475, 161 |
| 41 | Ginsenoside F3 | C41H70O13 | 815.4788 | -1.40 | 13.04 | m/z 161, 391, 475, 637, 769 |
| 42 | Ginsenoside F3 (isomer) | C41H70O13 | 815.4788 | -1.44 | 11.22 | m/z 161, 391, 475, 637, 769 |
| 43 | Pseudo-ginsenoside RT1 | C47H74O18 | 971.4846 | 0.40 | 9.64 | m/z 161, 763 |
| 44 | Ginsenoside Rs3 | C44H74O14 | 871.5050 | 0.40 | 25.25 | m/z 161, 459, 621,783 |
| 45 | Ginsenoside R2 | C41H70O13 | 815.4788 | -1.40 | 10.52 | m/z 161, 391, 475, 637, 769 |
| 46 | Ginsenoside Ra3 | C59H100O27 | 1285.6423 | -2.09 | 15.97 | m/z 1239, 1077, 945, 783, 621 |
| 47 | Ginsenoside Ra3 (isomer) | C59H100O27 | 1285.6423 | -2.10 | 17.52 | m/z 1239, 1077, 945, 783, 621 |
| 48 | Ginsenoside Rb3 | C53H90O22 | 1123.5895 | -1.00 | 18.45 | m/z 1077, m/z 1123 |
| 49 | Ginsenoside Rb3 (isomer) | C53H90O22 | 1123.5895 | -1.13 | 19.44 | m/z 1077, m/z 1123 |
| 50 | Ginsenoside Rk3 | C36H60O8 | 665.4259 | -1.30 | 20.68 | m/z 161, 619 |
| 51 | Protopanaxatriol | C30H52O4 | 521.3837 | -1.00 | 21.26 | m/z521,475,391 |
| 52 | Protopanaxatriol(isomer) | C30H52O4 | 521.3837 | -1.03 | 22.05 | m/z521,475,391 |
| 53 | Ginsenoside Ro | C48H76O19 | 1001.4952 | 2.40 | 25.07 | m/z 955,793,631,455 |

Conclusions

This study used the high-resolution SCIEX X500R QTOF System for identification of ginsenoside components. It uses SCIEX OS software along with the TCM MS/MS database for rapid, accurate identification of 53 ginsenoside components, showing strong resolving power and illustrating the benefits of the high-resolution SCIEX database in traditional Chinese medicine analysis. The high-resolution MS/MS TCM database contains almost a thousand TCM active ingredient MS/MS spectra; automatic data extraction can be used for matching and greatly decreases the identification time for Chinese medicines. It also allows for simple and accurate component identification.

The SCIEX X500R QTOF high-resolution system is the right tool when identifying Chinese medicine components. The IDA workflow can be used to ensure the integrity of the acquisition, and TOF MS and MS/MS data can be obtained for all components. The X500R's front end has all the advantages of

a triple quadrupole mass spectrometer, greatly improving its quantification capabilities, sensitivity, stability, and linear range.

SCIEX OS system software integrates instrument control, method editing, data acquisition, and reporting. It can perform simultaneous qualitative and quantitative analysis, wirelessly connect to other software, and simplify analytic workflow.

References

1. Shen Dan, Tang Shihuan, Lu Peng, et al. Ginseng-containing Chinese medicine pattern analysis [J]. Chinese Journal of Traditional Chinese Medicine, 2013, 38(12), 2027-2032.
2. Zhang Xiaoxu, Wang Hongping, Yang Yang, et al. Rapid identification of saponins in sunflower ginseng by ultra performance liquid chromatography tandem quadrupole time-of-flight mass spectrometry [J]. China Medical Journal, 2015, 12(9), 130-136.

AB Sciex is doing business as SCIEX.

© 2017 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-02-6090-A



Headquarters

500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales

For our office locations please call the division headquarters or refer to our website at sciex.com/offices